DE NasalOxy Nasal-/Oralbrille, 210 cm Schlauch
 EN NasalOxy nasal/oral cannula, 210 cm tubing
 FR NasalOxy lunettes naso-buccales, tubulure de 210 cm
 IT NasalOxy cannula oro-nasale, tubo da 210 cm
 ES NasalOxy cánula naso-oral, tubo de 210 cm
 PT NasalOxy cânula naso-oral, tubo de 210 cm
 NL NasalOxy neus-mondcanule, slang van 210 cm
 DK NasalOxy nasal/oral kanyle, 210 cm slange
 NO NasalOxy nasal/oral kanyle, 210 cm slange
 SE NasalOxy nasal/oral kanyl, 210 cm slang
 FI NasalOxy nenä-suukanyyli, 210 cm letku
 PL NasalOxy kaniula nosowo-ustna, rurka 210 cm
 RO NasalOxy canulă nazo-orală, tubulatură de 210 cm
 CN 鼻/口氧管, 210 厘米管子
 JP 鼻/口腔カニューレ、210cm チューブ
 KO 비강/구강 캐놀라, 210cm 튜브
 AR أنبوب طوله 210 سم NasalOxy، قنية أنفية فموية
 FA NasalOxy کانونل بینی / دهان، لوله 210 سانتی متر

Label Version: 02/2026



REF NC-290

LOT 2026004

 2026-03

 2031-02

(01)04262356150808
 (10)2026004
 (11)260301
 (17)310228

CNSAC MedShop GmbH
 Biebelrieder Str. 42
 97288 Theilheim
 Germany

Somnovum
 Industriestrasse 14
 CH-5036 Oberentfelden

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1 Operation

Look at the diagrams for how to put on the nasal/oral cannula:

1. Arrange the nasal/oral cannula as shown
2. Put the nasal/oral cannula in position
3. Secure under the chin

Proceed in reverse sequence to remove the nasal/oral cannula.

2 Intended use

The NasalOxy® nasal/oral cannula is intended for low-flow oxygen delivery (1-6 L/min) to spontaneously breathing patients weighing more than 30 kg. The device is designed for single patient use and facilitates oxygen delivery through both the nose and / or mouth. It is intended for use in clinical and homecare settings.

When used in clinical facilities without a change of patient, the nasal/oral cannula is replaced after 72 hours. In the domestic environment, the nasal/oral cannula is intended for a single patient with a maximum duration of use of 30 days.

2.1 Contraindication

The device is not intended for use:

- with oxygen supply systems using demand valves
- with high-flow oxygen therapy systems
- in applications requiring precisely controlled oxygen concentration (FiO₂)

2.2 Side effects

- mild skin irritation in the contact area between the upper lip and the base of the nose
- local discomfort due to prolonged use

3 Description of function

When oxygen is supplied and the nasal/oral cannula is attached in the area of the top lip, it generates a cloud of oxygen which can be inhaled through the nose and / or mouth. This method is scientifically proven; there is evidence of an adequate therapeutic effect even if the nasal passages are blocked.

4 Safety

In the EU: As a user and/or patient, you must report any serious incidents occurring in conjunction with the product to the manufacturer and to the responsible authority.

5 Product description

A diagram of the individual parts can be found on the title page.

1. Connection to oxygen supply
2. Attachment sleeve
3. Double tubing
4. Base plate with ball joint tube connections

6 Cleaning

It is not possible to clean or disinfect the entire nasal/oral cannula reliably. The manufacturer recommends cleaning just the outside of the nasal/oral cannula without its base plate under running water or by wiping over with a solution of alcohol-based disinfectant.

You can clean the base plate without ball joint tube connections as outlined below.

1. Wash your hands before cleaning.
2. Pull both ball joint tube connections off the base plate.
3. Wash the base plate in warm water and a mild detergent.
4. Wash the base plate by hand in warm water (approx. 30 °C) and with mild detergent (1 ml detergent to 1 l water) for 15 minutes.
5. Rinse the base plate with clean water.
6. Allow the base plate to air-dry.
7. Subject the base plate to a visual inspection.
8. Push both ball joint tube connections onto the base plate.

7 Disposal

In a domestic environment, you can dispose of all parts in domestic waste.
In a clinical environment, dispose of all parts in line with institutional regulations.

8 Technical Data

Oxygen flow range	1 - 6 L/min
Tubing length	210 cm (~ 7 ft)
Storage temperature	-20 °C bis +50 °C
Operating temperature	+5 °C bis +40 °C

9 Materials

Part of the nasal/oral cannula	Material
Patient circuit and connections	PVC (polyvinyl chloride)
Attachment sleeve, ball joints and base plate	PE (polyethylene)

10 Estimated Inspired Oxygen Concentration (FiO₂)

Estimated inspired oxygen concentration (FiO₂) depends on the patient's breathing pattern, tidal volume and inspiratory flow. Typical values during low-flow oxygen therapy may range approximately from 24% at 1 L/min to about 44% at 6 L/min. Actual oxygenation should be verified by monitoring the patient's oxygen saturation (SpO₂), e.g. by pulse oximetry.

O ₂ Flow (L/min)	Estimated FiO ₂ Nasal breathing	Estimated FiO ₂ Mouth breathing
1	~24–25 %	~23–24 %
2	~27–29 %	~26–28 %
3	~31–33 %	~30–32 %
4	~35–37 %	~34–36 %
5	~39–41 %	~38–40 %
6	~43–45 %	~42–44 %

Values are approximate and for guidance only.

11 Markings and symbols

The following symbols may be applied to the device, the device label, accessories or packaging..

Symbol	Description
	Unique Device Identifier (Unique Device Identification for medical devices)
	Permissible temperature range for transport and storage
	Use by the specified date
	Keep away from sunlight
	Order number
	Indicates that the product is a medical device
	Batch number

Symbol	Description
	Manufacturer
	Consult instructions for use
	Date of manufacture
	Single patient – multiple use
	CE marking (indicates conformity with applicable European directives/regulations)

**WARNING**

Risk of fire from the use of oxygen in combination with flammable substances. Supplying oxygen without a special safety device can lead to fire and injure people. Do not operate the nasal/oral cannula in the vicinity of excessive heat, naked flame or objects at risk of fire. Do not smoke or allow smoking near the nasal/oral cannula or oxygen supply. Oxygen-enriched environments significantly increase the risk of fire.

**WARNING**

Risk of infection if the nasal/oral cannula is reused! If the nasal/oral cannula is used for too long, infections may be transmitted. Regularly replace the used nasal/oral cannula with a new one before unhygienic contamination occurs.

**WARNING**

Risk of injury if the nasal/oral cannula malfunctions:
A damaged nasal/oral cannula or damaged individual parts may injure the patient.

- Only use undamaged nasal/oral cannulas.
- Do not crush or kink tubes, cannulas, and connections.
- Keep objects away on which the tubes may get caught.

**WARNING**

The effectiveness of oxygen therapy should be verified by monitoring oxygen saturation (SpO_2), for example with a pulse oximeter, especially when therapy is initiated or when flow settings are changed. Check all connections to ensure that there are no blockages or leaks in the device

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